## TROUBLE SHOOTING ZEPF NADELHALTER

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**24.950.10**

**24.950.12**

# TROUBLE SHOOTING ZEPF NEEDLE HOLDERS

## Maintenance and re-processing of needle holders

### BREAK OF JOINTS! CORROSION IN THE LOCK AND ON THE SURFACE!

Frequently needle holders are falsely claimed for the following reasons:

- Corrosion in the instruments joint, instrument breaks at the joint and at the tips.

In most cases (> 95%) we can exclude failures in material and workmanship. In almost all cases, the damage is caused by errors in the re-processing / hygiene chain. Below you will find instructions and explanations how to avoid these errors in the future.

### INSTRUMENT BREAK ON TIPS / JOINTS

Basically, all instruments need to be treated carefully. Do not throw them down or let them fall! Drops and mechanical stress due to impacts with other instruments may lead to hairline cracks. Particularly carbide inserts are very brittle and are very vulnerable for breaks by falls.

Please be sure you are using the right needle / thread size (type), which is recommended by the company for each kind of needle holder. E.g. micro needle holders are only recommended for sutures type 5/0 - 9/0.

All joint instruments need to be cleaned and sterilized in open position. This is particularly important for needle holders. When closed, the instrument tips are under contact pressure. Due to heating in the WD and during sterilization, the steel expands and the contact pressure on the instrument tip increases. This cannot occur in open position.

With Needle holders, especially for microsurgery with very fine points, a preparation in closed position can lead to stress cracking. In the beginning these cracks can absolutely be in the µm region, so that the final instruments break arises later, e.g. under stress.

Chlorides in varying concentrations in each drinking and mains water cannot be eliminated completely by a demineralization process of drinking water. While drying, chlorides focus punctually on the surface and cause pitting corrosion. To avoid this, the final rinse should be made with deionized water. These small corrosion holes, especially at the tips, inevitably lead to fractures.

### CORROSION ON JOINTS

If you ask about the meaning of terms such as stainless steel, it is often assumed that stainless steel is an indestructible, extremely durable material. The corrosion resistance of "STAINLESS STEEL" is based on the fact that they form passive layers on their surface, protecting the steel. Due to external influences, these protective passive layers can be destroyed, and this will lead to defects on the instruments.


# TROUBLE SHOOTING ZEPF NEEDLE HOLDERS

## Maintenance and re-processing of needle holders

### OILING AFTER CLEANING

After cleaning, the needle holders need to be cooled and lubricated before the function test to avoid metallic abrasion. Metal abrasion may cause the “seizure” of joints and lead to corrosion.

This is often a gradual process that finally causes a joint break of the needle holder. Joints of needle holders (scissors and forceps) must be oiled with a paraffin-based care product before the function test (e.g. REF **24.950.10** Medical White Oil).

#### Attention

Our **Double-Action Instruments** „**41.017.17** Micro Needle Holder“ and „**22.830.17** Micro Forceps“ need to be oiled additionally on both points of articulation (double transmission).

### REPROCESSING IN A WASHTRAY

We recommend a cleaning of our instruments in wash trays that ensure a fixation of the instruments and prevent from contact damage.

**MAKE SURE THAT YOUR VALUABLE INSTRUMENTS ARE ONLY USED FOR THEIR INTENDED APPLICATION!**

### RECOMMENDED CARE PRODUCTS

#### Cleaning Brushes with nylon bristles

**85.903.01** double-ended, sterilizable, packing unit: 3 pieces

**85.903.03** for aspirator Ø 1.5 mm

**85.903.04** for aspirator Ø 3 mm and Ø 5 mm

**Medical White Oil in a 50 ml dispensing bottle and in an Oil Pen, 12 ml**

**24.950.10**

**24.950.12**