## Mundspiegel

**Silberbeschichtung:** REF 24.061.22, 24.061.24

**Titanbeschichtung:** REF 24.062.24, 24.067.22, 24.067.24, 24.070.22

**Rhodiumbeschichtung:** REF 24.071.22, 24.071.24, 24.072.22, 24.072.24

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**Rhodium-coated:** REF 24.071.22, 24.071.24, 24.072.22, 24.072.24

### Dear customer!

You have opted for a premium product from our company – much appreciated.

For your continued satisfaction in the long-term we would like to go on to explain to you how to use this instrument:

### BASIC REMARKS

Thanks to the standardized thread-size, all the mirrors available from **HELMUT ZEPF** are compatible with all the handles supplied by **HELMUT ZEPF** and they can naturally all be disinfected and sterilized using customary procedures (autoclave max. 134 °C (273 °F) / hot-air 180 °C (356 °F)).

**Attention:** Do not use hard brushes, coarse sponges etc. in order to avoid scratching the glass surface or the coating.

Cleaning in the ultrasonic bath will not damage the mirror coating.

Acidic “refresh” products may only be used with rhodium mouth mirrors REF 24.071.22, 24.071.24, 24.072.22 and 24.072.24. As any silver coating can be destroyed by corrosive, acid “refresh” products, our rhodium-coated MEGA dental mirrors are appropriate if you use such products.

- In the ultrasonic or thermal disinfectant, water must not be contaminated by hard components or suspended solids. Hard particles (e.g. cement residue) in contaminated water can damage the mirror coating if you clean ultrasonically and clean in the thermal disinfectant.
- In the ultrasonic or thermal disinfectant other instruments must be stored in such a way to prevent them from touching the mirrors.
- Avoid short programs involving high and rapid temperature changes. Glass has a different coefficient of expansion to steel or plastic. A rapid and high change in temperature can induce tensions and result in a crack/fracture.

**Our rhodium mouth mirrors REF 24.071.22, 24.071.24, 24.072.22 and 24.072.24 are suitable for automated re-processing in mounted condition in our ZEPF **bionik** Universal Handles REF 26.194.01 - REF 26.194.10.**

### AVOID LIMESCALE AND ENCRUSTATION

1. Place instruments in a disinfection solution immediately after use. Then clean them in order to prevent any residual contamination.

2. Rinse with completely desalinated water. If you rinse with tap water, you have to dry the dental mirrors as otherwise salts and limescale can form on the instruments and water stains or white plaques can be deposited. These burn into the mirror surface at very high temperatures during sterilization and become permanent.

3. Descaling must function properly in thermal disinfectors. Otherwise the mirror will be obscured by a white, chalky plaque!

### REMOVING CHALKY WATER STAINS, MILKY WHITE LIMESCALE OR ENCRUSTATIONS

1. Remove with cleaning solution and a soft sponge cloth.

2. Remove with acidic products.

**This is only possible with rhodium-coated MEGA dental mirrors.**

This special mirror coating also acts as a non-stick layer like Teflon. As the rhodium coating is acid-proof, deposits can be removed by simply immersing in acidic refresh products (or acetic acid). Rhodium is an acid-proof precious metal like gold.

### Why isn't this possible for other mirrors?

No removal with acid products is possible for titanium- and silver-coated dental mirrors because lime burns into the glass surface at the high sterilization temperatures.

### INSTRUCTIONS FOR CLEANING AND STERILIZATION

#### Recommendation:

1. Decontamination immediately after use.
2. Cleaning: manual
3. Rinse with demineralized water and dry (meticulous inspection).
4. Packaging for sterilization.
5. Sterilization: Autoclave (max. at 134 °C / 273 °F), hot-air (hot-air 180 °C / 356 °F)

**Avoid using a thermal disinfectant to enable the mirror to have a long service life.**

#### Instructions:

Sterilization prior to initial use. Please follow the manufacturer's instructions for the device being used. Do not subject the instruments to a temperature of over 180 °C/356 °F (with hot-air sterilization). Only use instruments within their functional capacity.

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