# DENTAL-ZYLINDERAMPULLENSPRITZEN

**ZEPF REF 19.710.18**

## BESONDERE BEACHTUNG

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**ZEPF REF 19.710.18**

## SPECIAL OBSERVANCE

Please read the following information and application notes carefully. The product may only be used by experienced medical staff in accordance with these Instructions for Use. HELMUT ZEPF Medizintechnik GmbH does not give any recommendation regarding a method of treatment. The treating medical specialist staff is responsible for the manner of application and for the selection of the patient. If these instructions are not observed or if they are violated, this will lead to the expiry of the guarantee, and the safety of the patient will be endangered. If used in combination with further products, please always observe the indications for usage and the compatibility statements for such products. Please check the product and the packaging for completeness, intactness and the status of sterility before application. Do not use product if you have reason to doubt the completeness, the intactness and the non-hazardous status of sterility.

## PRODUCT SPECIFICATION / COMPATIBILITY

HELMUT ZEPF dental cartridge syringes are available in various designs, with swing-out mechanism to open the cartridge holder (the cartridge holder is firmly connected to the syringe), and also with slip-on device for separate cartridge holders (the cartridge holders can be replaced). The program includes models without aspiration, with self-aspiration, as well as with manual aspiration by means of hooklets or arrowhead.

## INDICATIONS

Infiltration- and nerve block anesthesia in dentistry

## CONTRAINDICATIONS

The contraindications to be observed are those of the applied anesthetic.

## COMPLICATIONS


No complications are known to occur under professional, proper application.

## WARNING INDICATIONS

Please observe when using models with aspiration hooklets! The aspiration hooklets **C** have to be in retracted position when inserting or removing the cartridges. They may otherwise be bent. Extend hooklets only to connect cartridge/plunger, please retract again when separating the connection! Retract the hooklets by turning the knurled screw **H** to the left hand side until the limiting stop becomes effective! Do not use the product in case of commone incompatibilities and / or interactions of the material.

- A** Screw thread for cannula.
- B** Cartridge holder with viewing windows.
- C** Three pivotable (swing-in and -out) aspiration hooklets at the end of the plunger provide for stable anchoring in the hollow stopper of the cartridge, and they ensure for smooth guidance when injecting and aspirating.
- D** Joint with swing-out mechanism (for opening the cartridge holder).
- E** Slip-on device for mounting the cartridge holder.
- F** Finger rest
- G** Plunger
- H** Knurled screw to extend and retract the hooklets. The hooklets are extended by turning the knurled screw **H** to the right-hand side (clockwise), and turning left (counter clockwise) will retract them. The screw-mechanism is equipped with a locking device. It is engaged as soon as the hooklets have been completely

extended by turning to the right-hand side. An unintentional separation of the connection is impossible. The effect of the locking device may be noticed as a slight initial resistance when turning the knurled screw to the left-hand side to separate the connection.

- J** Thumb plate ring; it permits the controlled application of pressure and it facilitates working "with delicate dexterity". The thumb plate ring makes a change of grip unnecessary, especially when pulling the plunger back for aspiration.
- K** Stationary aspiration arrowhead for solid stoppers.

## APPLICATION / SEQUENCE OF USE

- 1.** The cannula is tightly screwed onto the cartridge holder.
- 2.** Insert the cartridge
  - a)** For systems with swing-out mechanism: Pull the plunger (**G**) back completely. A spring resistance has to be overcome towards the end. Swing cartridge holder **B** out at a right angle. The holder is in open position. Insert the cartridge. Swing the holder back in alignment with the longitudinal axis. Thereby, completely retract the plunger backwards against the resistance of the spring. Release the plunger against the spring resistance. The cartridge will slip into the correct functioning position, centered in the cartridge holder by a guide ring.
  - b)** For systems with slip-on device: Pull the plunger (**G**) back completely. A spring resistance has to be overcome towards the end. Remove the cartridge holder **B** from the syringe by pulling it out of its rest laterally. Insert the cartridge into the cartridge holder. Slip the cartridge holder with the inserted cartridge into the rest of the syringe body. Thereby, completely retract the plunger backwards against the resistance of the spring. Release the plunger against the spring resistance. The cartridge will slip into the correct functioning position, centered in the cartridge holder by a guide ring.
- 3.** The plunger **G** can be advanced forward several millimeters against the inserted cartridge: The front end of the plunger with the hooklets **C** will be fitted tightly into the recessed void of the hollow stopper. The connection is created by turning the knurled screw (**H**) to the right-hand side (= extension of the hooklets). The locking device to safeguard the connection is engaged as soon as the hooklets have been completely extended. Alternatively: Anchor arrowhead in the rubber stopper of the cartridge by pushing the plunger forward.
- 4.** Aspiration: The aspiration device is essentially based on the creation of a stable connection between the hollow stopper of the cartridge and the plunger of the syringe, which will bear tractive force and force of pressure. This jam-free fastened connection permits the performance of aspiration and injection without problems. General complications relating to local anesthesia during dental treatment are often the result of an unintentional intravenous or intra-arterial injection. Failures of onset may also occur in case of an inadvertent intravascular injection, since the local anesthetic is carried away quickly from the location of the application. The aspiration check, performed "lege artis", will largely provide protection from injecting into a blood vessel.
- 5.** The removal of the cartridge is performed in reverse order.

## ASSEMBLY / DISMANTLING

see "APPLICATION"

# DENTAL CARTRIDGE SYRINGES

**ZEPF REF 19.710.18**

## CLEANING / STERILIZATION

After use, the dental cartridge syringe is to be sterilized through autoclaving before further use. The syringe is to be cleaned by suitable means before sterilization. If the free movement of the locking mechanism should perhaps be impaired despite regular cleaning, then it is recommendable to clean the syringe body in the ultrasonic cleaning bath.

### Cleaning

The methods, facilities and equipment validated by the user are always to be observed and checked for compatibility with the details given here in every case.

Reusable instruments and equipment must be checked regarding functionality and intactness immediately before and after use, as well as after cleaning operations.

Contaminated instruments are to be processed as soon as possible after use; dismantle the dismountable instrument into its components according to the dismantling instructions.

Rinse the instrument with clear, flowing water before machine-cleaning and disinfection. Remove coarse contaminations carefully by means of a soft cloth or soft brush. Clean manually or mechanically. Strictly avoid unwashed areas/product residue patterns. Always observe the manufacturer's instructions.

### Cleaning Process

Remove adhering contaminations using a soft brush with plastic bristles. Do not use scouring/abrasive cleaning agents or metal brushes. Clean lumen and passage ducts using soft, rounded brushes with plastic bristles. The lumen diameter and the diameter of the brush must correspond matchingly. Finally, rinse with distilled or demineralized, flowing water.

Dry instrument with absorbent, soft and lint-free cloth. Dry lumen and passage ducts using grease- and dust-free compressed air. Perform cleaning / disinfection only in a suitable cleaning- and disinfection machine, applying a correspondingly appropriate cleaning / disinfection method. Connect components with lumen and passage ducts directly to the special cleaning element of the machine.

To clean the instruments which have been provided with a cleaning duct and which cannot be dismantled, connect the instruments to the special cleaning element of the machine by means of Luer-Lock connector, if available. When selecting the washing program, take the material of the instrument to be cleaned into account (e.g. stainless instrument steel, chrome-plated surfaces, aluminum).

Follow the instructions of the appliance manufacturer. Finally, rinse with demineralized water. Allow for a sufficient drying phase.

Remove the instrument from the machine immediately after completion of the cleaning program. Clean thermo-resistant instruments using Vario-TD- or BGA-program.

Please only use alkaline cleaning agents. Do not use acidic cleaning agents.

### Sterilizing with Steam

Steam sterilization must be performed according to a validated steam sterilization process (e.g. sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN 17665-1).

If sterilization is to be performed using the fractionated vacuum method, then the 134° C/3 bar-program (273° F) is to be used, employing a minimum hold time of 4 minutes, drying time: 20 min. Allow equipment / instruments to cool down to room temperature before renewed use. After steam sterilization, store exclusively in suitable containers designated for such purpose.

## MAINTENANCE

Allow instruments to cool down to room temperature. Slightly lubricate movable parts with sterilizable, steam permeable servicing oil. Assemble dismountable instrument. Check instrument after every cleaning and disinfection with regard to: cleanness, function and intactness e.g. bent, broken, cracked, worn and broken off parts. Sort out and replace damaged or faulty instrument. If necessary, repeat the reprocessing. Use consumable parts only once. Check consumable parts for damage before every use.

HELMUT ZEPF recommends to handle the instruments prudently and carefully, as well as the unconditional observance of these directions for use in order to achieve greatest possible serviceable life. The serviceable life is influenced in a considerable degree by careful treatment and the application of adequate maintenance and cleaning measures.

## LIMITATION OF WARRANTY / DISCLAIMER

HELMUT ZEPF Medizintechnik GmbH guarantees to manufacture its products with greatest possible care. THIS IS THE ONLY VALID GUARANTEE, AND IT SHALL REPLACE ALL OTHER WARRANTIES GIVEN AND REPRESENTATIONS MADE. It shall be observed, that due to the biological differences of the persons to be treated, no product is always absolutely effective under all environmental conditions and circumstances. Components of the sets manufactured by HELMUT ZEPF Medizintechnik GmbH, as well as their individual components are compatible with each other. Before the use of individual products/sets of HELMUT ZEPF Medizintechnik GmbH in connection with products from other companies, the user must ensure the application-specific compatibility of the individual products.

HELMUT ZEPF Medizintechnik GmbH has no influence on the application of the product, on the diagnosis of the patient and on the handling of the product outside of the company. HELMUT ZEPF Medizintechnik GmbH can neither guarantee a beneficial nor a complication-free application of the product. HELMUT ZEPF Medizintechnik GmbH therefore assumes no liability for damages and costs. HELMUT ZEPF Medizintechnik GmbH will replace products showing a deficiency, which is to be represented by HELMUT ZEPF Medizintechnik GmbH. Employees of HELMUT ZEPF Medizintechnik GmbH shall not be authorized to amend the aforementioned conditions, to extend liability, or to accept or consent to additional product-related obligations.

## GENERAL NOTES

In accordance with U.S. law, this product may only be sold by a physician or upon his prescription.

All rights to change or modification of the product shall remain reserved. Products are free of latex. Made in Germany – manufactured by HELMUT ZEPF Medizintechnik GmbH.

## STORAGE

Please store the instruments in a dry, clean, and dust-free environment. Avoid large temperature fluctuations. This prevents the formation of condensate and related corrosion. In case of proper storage, there is no shelf-life specification or functional restriction for instruments supplied non-sterile after production.

## TRANSPORTATION

Transport in closed containers only.

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La HELMUT ZEPF Medizintechnik GmbH garantisce di fabbricare i propri prodotti con la massima cura possibile.