Skelettimplantate (Knochenstifte, Pins)

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Skeletal implants (bone pins, pins)

REF 47.560.03 - 47.561.06 Class: IIb R8 / UMDNS: 16-085

and accessories, instruments for insertion / use

REF 47.520.01

Dear valued customer!

You have decided for a high quality product from our company – many thanks.
In order to ensure your long-term satisfaction, we would like to explain the use of these implants and instruments as follows:

MANUFACTURER

HELMUT ZEPF Medizintechnik GmbH
Obere Hauptstr. 16-22
78606 Seitingen-Oberflacht, Germany
Telephone: +49 (0) 7464 / 98 88 0
Fax: +49 (0) 7464 / 98 88 88
E-Mail: info@zepf-dental.com
Internet: www.zepf-dental.com

Attention

Please read the information in this user information thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users. Please also note the inscriptions on the label.

INTENDED PURPOSE / INDICATION

The intended purpose of the products is rapid and complete restoration of the damaged bone function. Osteosynthesis implants cannot replace a broken bone, but can temporarily take over the mechanical function by keeping the fracture fragments in the anatomically correct position. The implant is tapped into the bone to fix inserted meshes or membranes, to isolate the bone from the soft tissue during regeneration and/or to fix the augmented bone mash which consists of autologous bone material and/or bone mixture or bone graft substitute. The pin implants can be combined with compatible products.

The implants are made of a TiAl6V4 alloy (ASTM F 136 in accordance with ISO 5832-3) and are biocompatible, as well as CT and MRI compatible. Investigation of CT and MRI compatibility has not been performed, however.

Contraindication

An insufficient amount or quality of bone mass to allow stable fixation. Infection or ulcers where the implant is to be placed.

Attention

Instruments and especially implants may only be used by persons who are specially trained or instructed for this purpose. The pins have to be applied with a compatible instrument. We recommend a pin applicator from Helmut Zepf Medizintechnik GmbH, art. no.: 47.520.01. It is recommended to explant the implants after the patient's recovery. Upon evaluation by the physician, the physician can decide under his sole responsibility whether the implants can remain in the patient's body.

APPLICATION & SAFETY INFORMATION

It is very important to inspect every product prior to use for visible damage and wear, e.g. cracks, breaks or defects. Never use damaged products.

SURFACE OF ALL IMPLANTS

The surface of all implants is chemically passive; they are not magnetic.

CORRECT SELECTION OF THE IMPLANTS

We refer to the fact that implants can only fulfill their function if the following basic rules are observed:

- When selecting the implants / prosthetics and fixators, care must be taken that they are selected according to the weight and level of activity of the patient and the bone fracture to be treated.
- Check that a suitable choice of biomechanics is made such that the forces transmitted by the implants are kept low.
- Severe deformations of the implants have to be avoided.
- Multiple deformations should always be avoided.
-  The reuse of implants is not permitted.

- It is absolutely necessary to inform the patient of the advantages and disadvantages of implants!

Due to the limited stability of the implants, avoid excessive strain from the patient's body weight. Failure to observe this precaution can result in severe consequences.

RECONDITIONING

The implants themselves may only be used once under any circumstances. As soon as an application is concluded, the implant must not be used again.

Following treatment of a patient infected with Creutzfeldt Jakob disease, the instruments and accessories for insertion / use must not be reused or reconditioned either. The instruments have to be disposed of.

GENERAL PRINCIPLES

All products must be cleaned, disinfected and sterilized before each use; this particularly applies to the first use after delivery, as all products are delivered unsterilized (cleaning and disinfection after removing the transport protection packaging; sterilization after packaging). A thorough cleaning and disinfection is an indispensable requirement for effective sterilization.

Skeletal implants (bone pins, pins)

Please note as part of your responsibility for the sterility of the products during use that

- generally, only suitable equipment and product-specific validated procedures are to be used for cleaning/disinfection and sterilization,
- the equipment used (WD, sterilizer, etc.) are to be regularly maintained and inspected, and
- the validated parameters are to be observed for each cycle.

Please also adhere to the legal requirements applicable in your country as well as the hygienic requirements of the medical practice or hospital. This particularly applies to the various requirements (e.g. in Germany according to attachment 7 of KRINKO RKI BfArM recommendation for processing) regarding to an effective prion inactivation (not applicable for the USA).

Remark:

Processing must be performed only by qualified staff in the central sterilization service department of the hospital or in the processing room of the medical practice. Hospital or medical practice are responsible for selection and application of required of protective equipment and hygienic measures.

PREPARATION AND TRANSPORT

The implants must be stored and transported to the preparation site in a closed container to avoid damage to the implants and/or the packaging and contamination of the environment.

CLEANING AND DISINFECTION

Principles

For cleaning and disinfection, if possible an automated procedure [WD (washer-disinfector)] should be used. A manual procedure – even using an ultrasound bath – should only be used according to country specific requirements (e.g. in Germany for critical B products automated procedure binding) and if an automated procedure is not available due to the significantly lower effectiveness and reproducibility.

AUTOMATED CLEANING AND DISINFECTION [WD (CLEANING AND DISINFECTION DEVICE)]

When selecting the WD, ensure that

- the WD generally has verified effectiveness (e.g. DGHM or FDA approval/clearance/registration or CE marking in accordance with DIN EN ISO 15883),
- if possible, a tested program for thermal disinfection (A0 value ≥ 3000 or – for older devices – at least 5 min at 90 °C/194 °F) is used (in chemical disinfection danger of disinfecting agent residues on the products),
- the program used is suitable for the products and contains sufficient rinsing steps (at least three degrading steps after cleaning (respectively neutralization, if applied) or conductance based rinsing control recommended in order to prevent effectively remnants of the detergents),
- for rinsing only sterile (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)¹ is used,

- air used for drying is filtered (oil-free, low-bacteria and low-particle) and
- the WD is regularly maintained, inspected, and calibrated.

When selecting the cleaning system, ensure that

- it is generally suitable for cleaning medical instruments made of metals,
- providing no thermal disinfection is used – a suitable disinfecting agent with verified effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is also used and that it is compatible with the cleaning agent used, and
- the chemicals used are compatible with the products (see “Material Stability” section).

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning agent and, if applicable, the disinfecting agent as well as specifications for rinsing must be adhered to.

Procedure:

1. Place the products into the WD. Ensure that the products do not touch.
2. Start the program.
3. Remove the products after the program has completed.
4. Inspect and pack the products as soon as possible after removal (see “Inspection” and “Packaging” chapters, possibly after additional drying in a clean area).

The verification of products’ general suitability for effective automated cleaning and disinfecting was provided by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory using the G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the Neodisher MediClean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). Here, the procedure described above was taken into consideration.

MANUAL CLEANING AND DISINFECTION

When selecting the cleaning and disinfecting agent, ensure that

- it is generally suitable for cleaning and disinfecting medical instruments made of metals,
- the cleaning agent is suitable for ultrasound cleaning (no foam formation),
- a disinfecting agent with verified effectiveness (e.g. VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this is compatible with the cleaning agent used, and
- the chemicals used are compatible with the products (see “Material resistance” chapter).

Combined cleaning/disinfecting agents should not be used if possible.

¹ In case of consideration of a lower water quality as sufficient based on the background of national recommendations (e.g. in Germany KRINKO RKI BfArM recommendation for processing).

Skeletal implants (bone pins, pins)

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning and disinfecting agent as well as specifications for rinsing, must be adhered to. Only use freshly prepared solutions, sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)¹ or only with a soft, clean, lint-free cloth (Attention: caution in case of products with rough surfaces, sharp edges or comparable aspects with danger of attachment of particles from the cloth!) and/or filtered air to dry.

Procedure: Cleaning

- Place the products for the predefined soaking time in the cleaning bath (in an ultrasound bath that is not already activated), so that the products are completely submerged. Ensure that the products do not touch. Support the cleaning by completely brushing all internal and external surfaces with a soft plastic brush.
- Activate the ultrasound for an additional minimum exposure time (but not less than 5 min).
- Then remove the products from the pre-cleaning bath and rinse them at least three times thoroughly (for at least 1 minute) with water.
- Inspect the products (see "Inspection" chapters).

Disinfection

- Place the inspected products in the disinfection bath for the predefined soaking time so that the products are completely submerged. Ensure that the products do not touch.
- Then remove the products from the disinfection bath and rinse them at least five times thoroughly (for at least 1 minute) with water.
- Dry the products with filtered compressed air.
- Pack the products as soon as possible after removal (see "Packaging" section, possibly after additional drying in a clean area).

The proof of the general suitability of the products for effective manual cleaning and disinfecting was provided by an independent, governmentally accredited and respected (§ 15 (5) German Law for Medical Devices) test laboratory using the Cidezyme/Enzol pre-cleaning and cleaning agent and the Cidex OPA disinfecting agent (Johnson & Johnson GmbH, Norderstedt). Here, the procedure described above was taken into consideration.

Inspection

Check all products after cleaning or cleaning/disinfecting for corrosion, damaged surfaces, chippings, contaminants and stains as well as remove damaged products. Any products that are still contaminated must be cleaned again and disinfected.

Packaging

Sort the cleaned and disinfected products into the corresponding sterilization tray.

Please pack the products or the sterilization trays in sterilization containers in accordance with the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature stability up to at least 138 °C (280 °F) sufficient steam permeability)
- sufficient to protect the products or sterilization packaging from mechanical damage
- undergo regular maintenance according to the manufacturer's specifications (sterilization containers)
- do not exceed a maximum weight of 10 kg per package/contents of the sterilization container.

STERILIZATION

For sterilization, only the following sterilization methods may be used; other sterilization methods are not allowed.

Steam sterilization

- Fractionated vacuum procedure^{2,3} (with sufficient product drying⁴)
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance assessment (PQ))
- maximum sterilization temperature 134 °C ((273 °F); plus tolerance in accordance with DIN EN ISO 17665)
- Sterilization time (exposure time at sterilization temperature):

Country	fractionated vacuum procedure
Germany	at least 5 min ⁵ at 134 °C (273 °F)
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min ³
France	at least 5 min ⁵ at 134 °C (273 °F) if required for prion inactivation, sterilization time 18 min
other countries	at least 5 min ⁵ at 132 °C (270 °F) / 134 °C (273 °F)

Gravity displacement is not recommended³.

Verification of the general suitability of the products for effective steam sterilization was provided by an independent, governmentally accredited and respected (§ 15 (5) MPG) test laboratory using the HST 6x6x6 steam sterilizer (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum procedure, as well as the instrument oil LAWTON MEDOIL. Here, the typical conditions in the clinic and medical practice and the procedure described above were taken into consideration.

The flash sterilization procedure is generally not permitted.

Do not use dry heat sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

² at least three vacuum steps

³ The use of the less effective gravity displacement is only permitted if the fractionated vacuum procedure is not available. It requires significantly longer sterilization times and must be validated by the user for each specific product, device, procedure and parameter.

⁴ The actual required drying time depends directly on the parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer status) and must, therefore, be determined by the user. Nonetheless, drying time should not be under 20 min.

⁵ or 18 min (prion inactivation, not relevant for the USA)

Skeletal implants (bone pins, pins)

Material Stability

When selecting the cleaning and disinfection agents, please ensure that they do not contain the following components:

- Organic, mineral and oxidizing acids (minimum permitted pH value 5.5)
- Alkalis/strong alkalis (neutral/enzymatic (max. permitted pH 8.5, mandatory requirement for products made of aluminum or other alkali-sensitive materials, see “Special Instructions” section) or alkaline cleaner (max. permitted pH 11, mandatory requirement for products with intended application in prion-critical areas, e.g. in accordance with Annex 7 of KRINKO RKI BfArM recommendation for treatment) recommended)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizing agents (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

Never clean products, sterilization trays or sterilization containers with metal brushes or steel wool.

All products, sterilization trays and sterilization containers can only be exposed to temperatures under 138 °C (280 °F).

STORAGE

Please store the instruments in a dry, clean, and dust-free environment. Avoid large temperature fluctuations. This prevents the formation of condensate and related corrosion.

In case of proper storage, there is no shelf-life specification or functional restriction for instruments supplied non-sterile after production.

After sterilization, the products must be stored dry and free of dust in the sterilization packaging.

REPAIRS

Repairs are not permissible under any circumstances.

INFORMATION ON VALIDATION OF RECONDITIONING

It is the responsibility of the user to validate their procedure appropriately.

HANDLING

Always treat all instruments and implants with the greatest care in transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

WARRANTY

Our company exclusively delivers inspected and flawless products to its customers. All our products are designed and manufactured such that the highest quality standards are fulfilled. Liability for products that are modified as compared with the original, misused or improperly treated or used, is excluded.

HELMUT ZEPF MEDIZINTECHNIK GMBH ACCEPTS NO LIABILITY IF THIS CUSTOMER INFORMATION HAS BEEN VERIFIABLY VIOLATED.

THESE PRODUCTS ARE ONLY INTENDED FOR SINGLE USE.

DECONTAMINATION INFORMATION

On the basis of statutory regulations and above all to protect our employees, we need a signed “Evidence of decontamination” with every return.

Please ensure that this “Evidence of decontamination” is enclosed filled out and signed with **every product return** (complaint / repair / other reason for return) and package the product such that there is no risk of injury for our personnel in goods receiving during unpacking.

This form is available for download on our website “www.zepf-dental.com”.

RELEASE VERSION – SEE FOOTER