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zepf-dental.com

AESTHETIC IS THE RESULT

CE 0297

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Schrauben

REF 47.572.50 - 47.572.55, 47.573.04S - 47.575.43

Klasse: IIb R8 / UMDNS: 16-101

GB  Application Guidance 7-12

Screws


FR  Instructions d'utilisation 13-18

Vis

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Tornillos

-  ▶ Ø 1,0 mm STAHL
-  ▷ Ø 1,2 mm CMS Titan
-  ▶ Ø 1,5 mm MICRO Titan
-  ▷ Ø 1,2 mm SCHIRM Stahl
-  ▷ Ø 1,4/1,9 mm EMERGENCY Titan

CMF SYSTEM
CRANIO MAXILLOFACIAL FRACTURE SYSTEM

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Schrauben

REF 47.572.50 - 47.572.55, 47.573.04S - 47.575.43 / Klasse: IIb R8 / UMDNS: 16-101

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
Dear customer,

You have decided for a high-quality product from our company – many thanks.
In order to ensure your long-term satisfaction, we would like to explain the use of these instruments as follows:

MANUFACTURER

HELMUT ZEPF Medizintechnik GmbH
Obere Hauptstr. 16-22
78606 Seitingen-Oberflacht, Germany
Phone: +49 (0) 7464 / 9888 0
Fax: +49 (0) 7464 / 9888 88
E-mail: info@zepf-dental.com
Internet: www.zepf-dental.com

Caution

Please read the information in this leaflet thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users. Please also note the inscriptions on the packaging.

INTENDED PURPOSE / INDICATION / CONTRAINDICATION

Intended purpose titanium screw / steel screw

The intended purpose is osteosynthetic treatment of ridge defects in the dental field. The titanium screws / steel screws are screwed into the bone to fix inserted membranes, titanium meshes or compatible bone blocks to isolate the bone during regeneration of the gingival soft tissue and/or to fix the augmented bone / bone replacement material taken from autologous/ allogeneic bone material or a mixture. For both materials, screwing into cortical bone and pre-drilling using compatible core drills is specified (see "CORRECT SELECTION OF SCREWS"). The products are used in maxillo-facial surgery. After healing, the products are explanted.

Contraindication

An insufficient amount or quality of bone mass to allow stable fixation. Infection or ulcers where the screws are to be placed. Material hypersensitivity, i.e., patient reactions to foreign bodies.

Material titanium

Titanium screws are made of a TiAl6V4 alloy (ASTM F 136 in accordance with DIN EN ISO 5832-3) and are biocompatible, as well as CT and MRI compatible.

Material implant steel

Steel implants are made of the implant material 1.4441 in accordance with DIN ISO 5832-1; they can generate artifacts in CT or MRI tomography.

To identify the materials, there are color indicators on the packaging of the screws. For more information, see "Color coding of the screws".

Undesirable side effects

Patients with a metal incompatibility (titanium or implant steel) may be irritated by contact with the product due to hypersensitivity. If these side effects occur, the procedure must be discontinued and if necessary, the required steps and measures must be taken.

Caution

Instruments and especially screws may only be used by persons who are specially trained or instructed for this purpose. Pre-drilling using compatible core drills is specified for screwing the screws into cortical bone.

In the case of steel implants, nuclear spin tomography can cause the material to heat up. This occurs very rarely, however, and the necessity of performing an MRI is the responsibility of the user.

APPLICATION & SAFETY INFORMATION

It is very important to inspect every product prior to use for damage and wear, e.g., cracks, breaks or defects. Never use damaged products.

SURFACE OF ALL SCREWS

The surface of the screws is chemically passive and not magnetic.

CORRECT SELECTION OF THE SCREWS

Please note that screws can only fulfill their function if the following basic rules are observed:

- When selecting the screws, care must be taken that they are selected according to the bone defects to be treated.

- The reuse of screws is not permitted.

- It is absolutely necessary to inform the patient of the advantages and disadvantages of the products!

Color coding of the screws

The identification of the screws is simplified by a colored indicator. The indicator is delivered together with the screws in a transport box. The diameter of the screws is indicated by the color of the indicator. The following color codes correspond to the same diameter:

- green → Ø 1.0 mm implant steel
- silver → Ø 1.2 mm titanium/implant steel
- red → Ø 1.4 mm titanium
- blue → Ø 1.5 mm titanium
- red → Ø 1.9 mm titanium

The indicators are labeled with the specific length of the screws on the top of the indicator. This does not release the user from the need to carry out length measurement prior to implantation. For each screw diameter and screws length there is an emergency screw, depending on the implant diameter and implant length, which is used if the bone density is lower than expected or depending on the user's indication. These products are marked with a red indicator.

Example indicator for "red":

Screws

Combination instruments and instruments for screw insertion

The color coding in the following table provides information about the combination with the corresponding screws.

Figure	REF and description
	<ul style="list-style-type: none"> 47.502.01 Screwdriver, small for Ø 1.0 mm, Ø 1.2 mm, Ø 1.4 mm, Ø 1.5 mm und Ø 1.9 mm
	Blades short <ul style="list-style-type: none"> 47.523.12* Ø 1.0/1.2/1.4 mm, 15 mm 47.523.15* Ø 1.5/1.9 mm, 15 mm
	Blades long <ul style="list-style-type: none"> 47.524.12* Ø 1.0/1.2/1.4 mm, 60 mm 47.524.15* Ø 1.5/1.9 mm, 60 mm
	Insertion aids <ul style="list-style-type: none"> 47.525.12* Ø 1.0/1.2/1.4 mm 47.525.15* Ø 1.5/1.9 mm
	<ul style="list-style-type: none"> 08.523.11* Gliding Hole Drill, Ø 1.1 mm, DD 9 mm 08.523.08* Core Hole Drill, Ø 0.8 mm, DD 7 mm
	<ul style="list-style-type: none"> 08.523.13* Gliding Hole Drill, Ø 1.3 mm, DD 12 mm 08.524.16* Gliding Hole Drill "Emergency", Ø 1.6 mm, DD 12 mm 08.523.09* Core Hole Drill, Ø 0.9 mm, DD 9 mm
	<ul style="list-style-type: none"> 08.524.16* Gliding Hole Drill, Ø 1.6 mm, DD 12 mm 08.524.19* Gliding Hole Drill "Emergency", Ø 2.0 mm, DD 12 mm 08.523.13* Core Hole Drill, Ø 1.3 mm, DD 12 mm

■ green → Ø 1.0 mm, □ silver → Ø 1.2 mm, ■ blue → Ø 1.5 mm,
■ red → Ø 1.4/1.9 mm, appropriate emergency screw
DD = Drilling Depth

Caution

A COMBINATION OF STEEL AND TITANIUM IS NOT PERMITTED!

The screws from Helmut Zepf Medizintechnik GmbH are supplied non-sterile. The products have to be firstly cleaned and sterilized prior to use.

RECONDITIONING

The screws themselves may only be used once under any circumstances. As soon as an application is concluded, the product must not be used again.

Following treatment of a patient infected with Creutzfeldt-Jakob disease, the instruments and accessories for insertion/use must not be reused or reconditioned either. The instruments have to be disposed of.

GENERAL PRINCIPLES

All products must be cleaned, disinfected and sterilized before each use; this applies in particular to the first use after delivery, as all products are delivered non-sterile (cleaning and disinfection after removal of the transport protective packaging; sterilization after packaging). Thorough cleaning and disinfection is an indispensable requirement for effective sterilization.

Please note as part of your responsibility for the sterility of the products in use that

- only suitable equipment and product-specific validated procedures are ever used for cleaning/disinfection and sterilization
- the equipment used (washer-disinfector, sterilizer, etc.) is to be regularly maintained, inspected and calibrated and
- the validated parameters are observed for each cycle.

Please also adhere to the legal requirements applicable in your country as well as the hygiene provisions of the medical practice or hospital.

This applies in particular to the different specifications (e.g., in Germany according to Annex 7 of the recommendation for reprocessing from KRINKO (German Commission for Hospital Hygiene and Infection Prevention), RKI (Robert Koch Institute), BfArM (Federal Institute for Drugs and Medical Devices)) regarding effective prion inactivation (not applicable for USA).

Note:

Reconditioning may only be carried out by trained staff in the central sterilization service department of the hospital or in the reconditioning room of the medical practice. The hospital or medical practice is responsible for the selection and application of the necessary protective equipment and hygiene measures.

CLEANING AND DISINFECTION

Principles

For cleaning and disinfection, if possible, an automated procedure (washer-disinfector or WD) should be used. A manual procedure – even using an ultrasonic bath – should only be used if an automated procedure is not available or in accordance with country-specific requirements (e.g., in Germany the automated procedure is mandatory for critical B-products) due to the significantly lower effectiveness and reproducibility.

AUTOMATED CLEANING / DISINFECTION

WD (washer-disinfector)

When selecting the WD ensure that

- the WD generally has verified effectiveness (e.g. The German Society for Hygiene and Microbiology (DGHM) or FDA approval/cleanance/registration or CE marking according to DIN EN ISO 15883)
- if possible, a tested program for thermal disinfection (A0 value > 3000 or – for older devices – at least 5 min. at 90°C/194°F) is used (in chemical disinfection, danger of disinfecting agent residues on the products)
- the program used is suitable for the products and includes sufficient rinsing steps (at least three degrading steps after cleaning (or neutralisation if applied) or conductance-based rinsing control recommended to effectively prevent detergent residues)
- only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water¹) is used
- the air used for drying is filtered (oil-free, low-bacteria and low-particle) and
- the WD is regularly maintained, checked and calibrated.

¹ Should you consider a lower water quality as sufficient against the background of national recommendations (e.g. KRINKO/RKI/BfArM recommendation for treatment), this is your sole responsibility.

Screws

When selecting the cleaning system to be used, ensure that

- it is generally suitable for cleaning invasive medical devices made of metals
- providing no thermal disinfection is used - a suitable disinfecting agent with verified effectiveness (e.g. The German Association for Applied Hygiene (VAH) / The German Society for Hygiene and Microbiology (DGHM) or FDA/EPA approval/cleanance/registration or CE marking) is also used and that this is compatible with the cleaning agent used and
- the chemicals used are compatible with the products (see chapter "Material Stability").

The concentrations, temperatures and soaking times as well as instructions for rinsing given by the manufacturer of the cleaning agent and, if applicable, disinfectant must be strictly adhered to.

Procedure:

1. Place the products in the WD. Ensure that the products do not touch each other.
2. Start the program.
3. Remove the products from the WD at the end of the program
4. Inspect and pack the products as soon as possible after removal (see chapters "Inspection" and "Packaging", possibly after additional drying in a clean area).

The verification of the products' general suitability of the products for effective automated cleaning and disinfection has been provided by an independent, officially accredited and recognized (§ 15 (5) German Medical Devices Act (MPG)) test laboratory using the WD G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the neodisher MediClean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). Here, the procedure previously described was taken into consideration.

MANUAL CLEANING AND DISINFECTION

When selecting the cleaning and disinfecting agent to be used, ensure that

- it is generally suitable for cleaning and disinfecting invasive medical devices made of metals
- the cleaning agent is suitable for ultrasonic cleaning (no foam formation)
- a suitable disinfecting agent with verified effectiveness (e.g. VAH/DGHM or FDA/EPA approval/cleanance/registration or CE marking) is also used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the products (see chapter "Material Stability").

Combined cleaning/disinfecting agents should not be used if possible.

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning and disinfecting agent, must be strictly adhered to. Only use freshly prepared solutions, only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)¹ and only a soft, clean and lint-free cloth and/or filtered air for drying.

Caution

The products have a rough surface and sharp edges on which particles from the cloth can be caught!

Procedure: Cleaning

1. Place the products in a sufficiently large pre-cleaning bath (in an ultrasonic bath that has not been activated yet) for the specified soaking time so that the products are completely submerged. Ensure that the products do not touch each other. Support the cleaning by completely brushing all internal and external surfaces with a soft plastic brush.
2. Activate the ultrasound for another minimum soaking time (but not less than 5 min.).
3. Then remove the products from the cleaning bath and rinse thoroughly (for at least 1 min.) with water at least three times.
4. Inspect the products (see chapter "Inspection").

Procedure: Disinfection

5. Place the products for the predefined soaking time in the cleaning bath so that the products are completely submerged. Ensure that the products do not touch each other.
6. Then remove the products from the disinfection bath and rinse thoroughly (for at least 1 min.) with water at least five times.
7. Dry the products by blowing off/out with filtered compressed air.
8. Pack the products as soon as possible after removal (see chapter "Packaging", possibly after additional drying in a clean area).

The verification of the products' general suitability for effective automated cleaning and disinfecting were provided by an independent, officially accredited and recognized (§ 15 (5) German Medical Devices Act (MPG)) test laboratory using the pre-cleaning and cleaning agent Cidezyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt). Here, the procedure described above was taken into consideration.

Inspection

Check all products after cleaning or cleaning/disinfection for corrosion, damaged surfaces, disinfection, contaminants as well as stains and discard damaged products. Products that are still contaminated must be cleaned again and disinfected.

Packaging

If necessary, sort the cleaned and disinfected products into the corresponding sterilization tray.

Please pack the products or the sterilization trays in sterilization containers in accordance with the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature stability up to at least 138°C (280°F) sufficient steam permeability)
- sufficient to protect the products or sterilization packaging from mechanical damage
- regular maintenance according to the manufacturer's specifications (sterilization container)
- Do not exceed a maximum weight of 10 kg per package/content of the sterilization container.

¹ Should you consider a lower water quality as sufficient against the background of national recommendations (e.g. KRINKO/RKI/BfArM recommendation for treatment), this is your sole responsibility.

Screws

STERILIZATION

For sterilization, only the following sterilization methods may be used; other sterilization methods are not permitted.

Steam sterilization

- Fractionated vacuum procedure^{2,3} (with sufficient product drying⁴)
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Validation in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance qualification (PQ))
- Maximum sterilization temperature 134°C (273°F); plus tolerance in accordance with DIN EN ISO/ANSI AAMI ISO 17665
- Sterilization time (soaking time at sterilization temperature):

Country	fractional vacuum procedure
Germany	at least 5 min. ⁵ at 134°C (273°F)
USA	min. 4 min. at 132°C (270°F), drying time min. 20 min. ³
France	min. 5 min. ⁵ at 134°C (273°F) if required for prion inactivation, sterilization time 18 min.
other countries	at least 5 min. ⁵ at 132°C (270°F) / 134°C (273°F)

Gravity displacement is not recommended³.

The verification of the products' general suitability of the products for effective steam sterilization has been provided by an independent, officially accredited and recognized (§ 15 (5) MPG) test laboratory using the HST 6x6x6 steam sterilizer (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum procedure and LAWTON MEDOIL (oils for joints and friction surfaces).

Here, typical conditions in clinics and medical practices and the procedure described above were taken into consideration.

The flash sterilization process is generally not permitted.

Do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

Material Stability

When selecting cleaning and disinfecting agents, please ensure that they do not contain the following components:

- Organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- Alkalis/strong alkaline solutions (neutral/enzymatic (max. permissible pH 8.5, mandatory requirement for products made of aluminum or other alkali-sensitive materials or alkaline cleaner (max. permitted pH 11, mandatory requirement for products with intended use in prion-critical areas, e.g. In accordance with Appendix 7 of the KRINKO RKI BfArM recommendation)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizing agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

Never clean products, sterilization trays or sterilization containers with metal brushes or steel wool.

All products, sterilization trays or sterilization containers can only be exposed to temperatures under 138°C (280°F).

STORAGE

Store screws in a dry, clean and dust-free environment. Avoid large temperature fluctuations. This prevents the formation of condensate and the related corrosion. A shelf-life specification or functional restriction is not given for screws supplied non-sterile after manufacture if they are stored properly.

After sterilization, the products must be stored dry and free of dust in the sterilization packaging.

REPAIRS

Repairs are not permissible under any circumstances.

INFORMATION ON VALIDATION OF RECONDITIONING

It is the responsibility of the user to validate their procedure appropriately.

HANDLING

Caution

Always treat all screws with the greatest care in transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

WARRANTY

Our company exclusively delivers inspected and flawless products to its customers. All our products are designed and manufactured such that the highest quality standards are fulfilled. Liability for products that are modified as compared with the original, misused or improperly treated or used, is excluded.

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ACCEPTS NO LIABILITY IF THIS CUSTOMER INFORMATION HAS BEEN VERIFIABLY VIOLATED.

THESE PRODUCTS ARE ONLY INTENDED FOR SINGLE USE.

DECONTAMINATION INFORMATION

Information on return deliveries

On the basis of statutory regulations and above all to protect our employees, we need a signed "Evidence of decontamination" with every return. Please ensure that this "Evidence of decontamination" is enclosed filled out and signed with every product return (complaint/repair/other reason for return) and package the product such that there is no risk of injury for our personnel in Goods Receiving during unpacking.

The form is available for download on our website at "www.zepf-dental.com".

Attention

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A reporting form for this purpose can be found on our website.

RELEASE VERSION – SEE FOOTER

² at least three vacuum steps


³ The use of the less effective gravity displacement is only permitted if the fractionated vacuum method is not available. It requires significantly longer sterilization times and must be validated by the user for each specific product, device, procedure and parameter.

⁴ The actual required drying time depends directly on parameters that are the sole responsibility of the user (loading configuration and density, sterilizer status, ...) and must, therefore, be determined by the user. Nevertheless, drying times should not be under 20 min.

⁵ or extended sterilization time (e.g. 18 min.) for prion inactivation according to national regulations (not relevant for USA)

Screws

Color Pins

- 1.0 mm = "green"
- 1.2 mm = "silver"
- 1.5 mm = "blue"
- 1.4 mm = "red"
- 1.9 mm = "red"

Color Pins are included in each packing unit of screws. An easy identification of the screws is guaranteed by using the length markers in the tray. The colored pin with length marking is assigned to the respective screw diameter, thus allowing the determination of the screw length and a clear lot traceability in the box.


Visit our Product Assistant

www.zepf-dental.com/shop/en/cmff

Ø 1.0 mm	
STEEL Screws	PU: 5 pieces
47.573.04S*	Ø 1.0 x 4 mm
47.573.06S*	Ø 1.0 x 6 mm
47.573.08S*	Ø 1.0 x 8 mm
47.573.10S*	Ø 1.0 x 10 mm

Ø 1.2 mm	
CMS Titanium Screws	PU: 5 pieces
47.574.04*	Ø 1.2 x 4 mm
47.574.06*	Ø 1.2 x 6 mm
47.574.08*	Ø 1.2 x 8 mm
47.574.10*	Ø 1.2 x 10 mm
47.574.12*	Ø 1.2 x 12 mm

EMERGENCY Screws Ø 1.4 mm / PU: 3 pieces	
47.574.35*	Ø 1.4 x 5 mm
47.574.40*	Ø 1.4 x 10 mm


Ø 1.5 mm	
MICRO Titanium Screws	PU: 5 pieces
47.575.04*	Ø 1.5 x 4 mm
47.575.06*	Ø 1.5 x 6 mm
47.575.08*	Ø 1.5 x 8 mm
47.575.10*	Ø 1.5 x 10 mm

EMERGENCY Screws Ø 1.9 mm / PU: 3 pieces	
47.575.36*	Ø 1.9 x 6 mm
47.575.39*	Ø 1.9 x 9 mm
47.575.43*	Ø 1.9 x 13 mm

ZEPF Umbrella Screws for the CMF SYSTEM

Successful defect regeneration requires space and calm. Our umbrella screws create space and a stable area for the inserted augmentation material. The covering soft tissue is kept away and prevents the augmentation from "shifting/compressing" when particulate augmentation material is used.

The steel screws are available in different lengths with Ø 1.2 mm. The corresponding screwdriver blades are REF 47.524.12 or REF 47.523.12.

UMBRELLA Steel Screws Ø 1.2 mm / PU: 1 piece	
47.572.50*	Ø 1.2 x 8 mm, head Ø 4 mm
47.572.51*	Ø 1.2 x 10 mm, head Ø 4 mm
47.572.52*	Ø 1.2 x 12 mm, head Ø 4 mm
47.572.53*	Ø 1.2 x 8 mm, head Ø 6 mm
47.572.54*	Ø 1.2 x 10 mm, head Ø 6 mm
47.572.55*	Ø 1.2 x 12 mm, head Ø 6 mm

Explanation of symbols

	Consult instructions for use
	Manufacturer
	Date of manufacture
	Use-by-date
	LOT no. (batch no.)
	Not for reuse – for single use only
	Non-sterile

	Unique device identifier
	Quantity
	HIBC Barcode
	CE marking
	Keep dry
	Caution
	Medical device

0703091032_47_SCHRAUBEN_HZ_BA | 01.02.2024 | 12/24

78606 Seitingen-Oberflacht, Germany | Tel.: +49 (0) 74 64 / 98 88 -0 | Fax: +49 (0) 74 64 / 98 88 -88 | info@zepf-dental.com | www.zepf-dental.com


Medical devices indicated here bear the CE mark. * Product higher than class I. Notified body will be indicated in the delivery note.

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